G I B A U D

BANDAGE HERNIAIRE HERNIA BELT

6089

BANDAGE HERNIAIRE FORT STRONG HERNIA BELT

6351

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EN | Instructions for Use

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SK | Návod na použitie

SL | Navodila za uporabo

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
If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

INDICATION

Uni- or bilateral inguinal reducible hernia, pre- and/or post-operative wear.

PERFORMANCE

The device provides localized pressure. The device reduces hernia symptoms during use.

CONTRAINDICATION

Incarcerated hernia (non-reducible).

WARNINGS • PRECAUTIONS OF USE

This product contains Natural Rubber Latex which may cause allergic reactions. When adjusting the strap, tighten it so that it is firm but comfortable.

Do not wear for extended periods without medical follow-up.

The device must not be used in direct contact with a damaged skin.

The device is a single patient use device. Do not use for multiple patients.

If any problems occur with the use of the device, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional.

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

FEATURES OF THE DEVICE

Anatomical pads (10.5 x 6 cm).

Removable leg straps for unilateral use.

Plastic buckles.

SET UP AND OPERATION

The first application should be carried out by a qualified healthcare professional to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend.

While lying down, position the pad and the belt over the hernia orifice (reduced hernias). Attach the belt and the leg straps, adjust the strap so that it is tight but comfortable. Check that the device is effective when standing; the hernia should be held by the pad.

If the hernia is unilateral, only use the pad on the side affected by the hernia. In this case, we recommend using both leg straps.

CARE INSTRUCTION, MAINTENANCE AND DISPOSAL

Keep the device away from sources of excessive heat.

Hand wash in warm soapy water (maximum 40 ° C), rinse thoroughly.


Wring out gently; do not bleach; dry flat, away from a source of excessive heat (ex: radiators), do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

MEASUREMENT AND SIZE CHART

The size is to choose in accordance with the waist circumference, measured 2 cm above the iliac crests (in cm).

1	2	3
76 - 89	89 - 104	104 - 117

If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

INDICATION

Uni- or bilateral inguinal reducible hernia, permanent wear.

PERFORMANCE

The device provides localized pressure. The device reduces hernia symptoms during use.

CONTRAINDICATION

Incarcerated hernia (non-reducible).

WARNINGS • PRECAUTIONS OF USE

This product contains Natural Rubber Latex which may cause allergic reactions. When adjusting the strap, tighten it so that it is firm but comfortable.

Do not wear for extended periods without medical follow-up.

The device must not be used in direct contact with a damaged skin.

The device is a single patient use device.

Do not use for multiple patients.

If any problems occur with the use of the device, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional.

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

FEATURES OF THE DEVICE

Anatomical pads (11 x 7 cm).

Removable leg straps for unilateral use.

Metal buckles.

SET UP AND OPERATION

The first application should be carried out by a qualified healthcare professional to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend.

While lying down, position the pad and the belt over the hernia orifice (reduced hernias). Attach the belt and the leg straps, adjust the strap so that it is tight but comfortable. Check that the device is effective when standing; the hernia should be held by the pad.

If the hernia is unilateral, only use the pad on the side affected by the hernia. In this case, we recommend using both leg straps.

CARE INSTRUCTION, MAINTENANCE AND DISPOSAL

Keep the device away from sources of excessive heat.


Hand wash in warm soapy water (maximum 40°C), rinse thoroughly. Wring out gently; do not bleach; dry flat, away from a source of excessive heat (ex: radiators), do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

MEASUREMENT AND SIZE CHART

The size is to choose in accordance with the waist circumference, measured 2 cm above the iliac crests (in cm).

1	2	3
75-90	90-105	105-120

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GIBAUD S.A.S.U