G I B A U D

LOMBOGIB® DOUBLE ACTION

6567/6658/6667/6666

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EN | Instructions for Use

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SV | Bruksanvisning

NO | Bruksanvisning


DA | Brugsanvisning

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
MISE EN PLACE
SET UP AND OPERATION
COME SI INDOSSA
ANLEGEN
PLACERING
PLASSERING
PÅSÆTNING
PAIKALLEEN ASETTAMINEN
ZAKŁADANIE


1

2

3

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
If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

INDICATIONS

Lower back conditions that may benefit from gross motion restriction, such as:

- Muscular low back pain (lumbago).
- Sciatica.
- Lumbosacral strains.
- Chronic low back pain.

PERFORMANCE

The device provides lumbar spine gross motion restriction.

The device reduces stress on the low back in flexion and extension.

CONTRAINDICATION

The devices must not be used by individuals for whom compression is contraindicated.

WARNING • PRECAUTION OF USE

Patients suffering from pelvic floor weakness should consult with their medical professional before using the orthosis.

The orthosis must not be used in direct contact with a damaged skin

This product contains Naturel Rubber Latex which may cause allergic reactions.

Do not wear for extended periods without medical follow-up

The orthosis is a single patient use device. Do not use for multiple patients.

If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Gradual withdrawal possibly combined with functional physiotherapy.

FEATURES AND BENEFITS

The device is a lumbar stretchable belt featuring four V-shaped posterior stays.

Ergonomic finger loops facilitate installation.

Double tightening: the large double straps are used to adapt the tightening according to activity and pain.

SET UP AND OPERATION

The first trial should be carried out by a qualified expert to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user, a family member or a friend.

During testing, the fitter must shape the V-shaped ribs of the brace to adjust it to your form.

- Wear over an undergarment or light clothing,

1) Position the bottom part of the belt at the top of the gluteal cleft.

2) Center the dorsal ribs on either side of the spine. The Gibaud logo must be on the right side; the ribs form a V on either side of the lumbar column.

3) The belt must be closed by stretching the two sides symmetrically using the finger holes.

CARE INSTRUCTION, MAINTENANCE AND DISPOSAL


Store the brace away from sources of extreme heat (including radiators and tumble dryers). Hand wash in warm soapy water (maximum 40 ° C), rinse thoroughly. Wring out gently; do not bleach; dry flat, do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

MEASUREMENT AND SIZE CHART

The device is available in two heights: 21 cm (6666/6667) & 26 cm (6567/6658)
The size is to choose in accordance with the waist circumference, measured 2 cm above the iliac crests (in cm)

	0	1	2	3	4	5
	70-80 cm	80-90 cm	90-100 cm	100-110 cm	110-125 cm	125-140 cm

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GIBAUD S.A.S.U