GIBAUD

LOMBOGIB® EVOLUTION

6660/6661/6669/6670

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EN | Instructions for Use

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SV | Bruksanvisning

NO | Bruksanvisning

DA | Brugsanvisning

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If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopedic consultant for advice.

INDICATIONS

Lower back conditions that may benefit from gross motion restriction, such as:

- Muscular low back pain.
- Sciatica.
- Lumbosacral strains.
- Chronic low back pain.

PERFORMANCE

The device provides lumbar spine gross motion restriction.

The device reduces stress on the low back in flexion and extension.

CONTRAINDICATION

Must not be used by individuals for whom compression is contraindicated.

WARNING • PRECAUTION OF USE

Patients suffering from pelvic floor weakness should consult with their medical professional before using the orthosis.

The orthosis must not be used in direct contact with a damaged skin.

This product contains Naturel Rubber Latex which may cause allergic reactions.

Do not wear for extended periods without medical follow-up.

The orthosis is a single patient use device. Do not use for multiple patients.

If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional.

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Gradual withdrawal possibly combined with functional physiotherapy.


FEATURES AND BENEFITS

Lombogib® Evolution is an extensible brace featuring V-shaped posterior stays.

A soft back pad can be added/removed. Additional stabilization is provided by lumbosacral reinforcement resting on the iliac crest. The brace is made of breathable material and available in two different heights: 21 and 26 cm.

MEASUREMENT AND SIZE CHART

The device is available in two heights: 21 cm (6669/6670) & 26 cm (6660/6661)
The size is to choose in accordance with the waist circumference, measured 2 cm above the iliac crests (in cm)

	0	1	2	3	4	5
	70-80 cm	80-90 cm	90-100 cm	100-110 cm	110-125 cm	125-140 cm

SET UP AND OPERATION

The first trial should be carried out by a qualified expert to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend.

During testing the fitter must shape the V-shaped stays of the brace to adjust it to your form.

1) Before fitting: open the belt, lay it down flat, remove both finger holes and separate the two parts of the belt.

2) Position the belt by placing the strap with the Gibaud logo on the right side and center the opening in the back piece over the backbone. The lower section of the back piece must be situated at the top of the buttock crease.

3) Close the belt without tightening it.

4) Find the finger holes and place your thumbs in the holes, hands facing forward; pull the straps towards the front to tighten firmly but comfortably. Attach the finger holes to the front of the belt.

The tightness of the belt can be adjusted according to pain and activity.

To remove the belt:

- Unhook the finger holes, leaving them unattached.

- Open and remove the belt.

CARE INSTRUCTION, MAINTENANCE AND DISPOSAL

Store the brace away from sources of extreme heat (including radiators and tumble dryers). Hand wash in warm soapy water (maximum 40 ° C), rinse thoroughly. Wring out gently; do not bleach; dry flat, do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

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GIBAUD S.A.S.U