G I B A U D

MALLÉOGIB[®] ARTICULÉE

ARTICULATED MALLÉOGIB[®]

(6585)

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EN | Instructions for Use

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If you have any doubts regarding the applications and use of this medical device, do not hesitate to ask your doctor and / or fitter, pharmacist, orthotist, or orthopaedic consultant for advice.

INDICATIONS

For conditions that may benefit from ankle inversion/eversion restriction, such as :

- Acute mild to severe ankle sprains.
- Chronic ankle instability.
- Post-cast and post walker support.
- Ankle injury rehabilitation and/or post-surgical use.
- Prophylactic use.

PERFORMANCE

The device provides inversion and eversion restriction to the ankle. The device improves proprioception.

CONTRAINDICATION

None known.

WARNINGS • PRECAUTIONS OF USE

The brace should not be used in cases where the injury requires relative or absolute immobilization or surgical procedure.

It's forbidden to drive a vehicle whilst wearing the brace.

Inspect the device carefully before applying it.

When adjusting the straps pay attention to have the tightness firm but comfortable.

Use footwear when walking outdoors or during sportive activities.

Avoid walking with the orthosis without footwear on slippery surfaces, particularly if the surface is wet (bathroom).

Do not wear for extended periods without medical follow-up.

The orthosis must not be used in direct contact with a damaged skin.

The orthosis is a single patient use device. Do not use for multiple patients.

If any problems occur with the use of the orthosis, e.g. pain or appearance of local signs, remove the orthosis and contact your healthcare professional.

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

FEATURES OF THE DEVICE

Combination of a heel cup to hold the rearfoot and two rigid shells.

Option to lock the joint straight during the acute phase.

Resilient padding.


SET UP AND OPERATION

The first application must be carried out by a qualified healthcare professional to understand the correct fitting and how to put the device on properly. Subsequently, the product should be put on the patient in the same way by the user or a family member or a friend.

1- Open the two straps and spread apart the shells.

MEASUREMENT AND SIZE CHART

Size of the orthosis to choose according to the shoe size (UK).

1	2	3
2-4	5-8	9-11

- 2- Place the heel into the opening on the cup.
- 3- Apply the shells on each side of the foot so that the malleoli are covered.
- 4- Tighten the lower strap and then the upper strap.
- 5- With the brace on, slide the foot into a spacious, stable shoe. Do not use pumps or any shoes with high heels. If the shoe has a thick, removable inner sole, taking it out will make walking more comfortable.
- 6- Loosening may occur with use. Retighten the straps if this occurs.
- 7- If you feel local discomfort and/or tingling, release the straps and reposition the shells.
- 8- On medical advice, the joint on the brace can be temporarily locked using the screws provided in the box.

CARE INSTRUCTION, MAINTENANCE AND DISPOSAL

Keep the device away from sources of excessive heat.

Wipe the plastic parts with a clean, damp cloth. The removable lateral protection can be machine washed (delicate cycle, low temperature).

Hand wash in warm soapy water (maximum 40°C), rinse thoroughly. Wring out gently; do not bleach; dry flat, away from a source of excessive heat (ex: radiators), do not tumble dry; do not iron; do not dry clean.

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

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GIBAUD S.A.S.U